# Application News

## No. C132

LC/MS  
 Liquid Chromatography Mass Spectrometry

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Comprehensive Analysis of Primary and Secondary Metabolites in Citrus Fruits Using an Automated Method Changeover UHPLC System and LC/MS/MS System [LCMS-8050]

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Multiple Reaction Monitoring (MRM) Analysis of Primary and Secondary Metabolites

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Comparison of Flesh and Peel Constituents

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表 2 目标成分的标准曲线线性范围  
Calibration Curve Ranges of Target Compounds

Amino acids	CAL of range	Organic acids	CAL of range
Cystine	1-100	Tartaric acid	50-10000
Aspartic acid	5-100	2-Ketoglutaric acid	10-10000
Asparagine	5-100	Isocitric acid	50-10000
Serine	5-100	Malic acid	10-5000
4-Hydroxyproline	1-100	Lactic acid	50-10000
Glycine	5-100	Citric acid	50-10000
Lysine	1-100	Pyroglutamic acid	10-10000
Cysteine	50-100	Succinic acid	10-1000
Threonine	5-100	Fumaric acid	500-1000
Glutamic acid	1-100	Maleic acid	50-10000
Alanine	5-100		
Proline	1-100		
Ornithine	5-1000		
Glutamine	5-100		
Histidine	5-100		
Arginine	5-100		
GABA	5-100		
Valine	1-100		
Methionine	5-100		
Tyrosine	5-100		
Isoleucine	5-100		
Leucine	10-100		
Phenylalanine	1-100		
Tryptophan	5-500		
Sugars	CAL of range		
Rhamnose	50-1000		
Fructose	50-5000		
Glucose	50-5000		
Sucrose	100-5000		
Maltose, Lactose	200-2000		
Carotenoids	CAL of range	Flavonoids	CAL of range
Fucoxanthin	0.1-100	Naringin	10-1000
Violaxanthin	1-100	Quercetin	5-1000
Astaxanthin	0.5-100	Kaempferol	5-1000
Lutein	0.1-100	Nobiletin	0.1-1000
Zeaxanthin	0.5-100	Tangeretin	0.1-1000
Canthaxanthin	0.05-100		
$\beta$ -Cryptoxanthin	0.05-100		
Lycopene	50-100		
$\beta$ -carotene	0.1-100		


图 5 果实与果皮的主成分分析 (PCA)  
Principal Component Analysis (PCA) of Flesh and Peel

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Comparison of Sugars (Flesh and Peel)

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
图6 果实与果皮所含糖的量与组成的比较  
Comparison of Sugar Quantity and Sugar Composition in Flesh and Peel

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Comparison of Secondary Metabolites (Flesh and Peel)

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
图8 果实的主成分分析 (PCA)  
Principal Component Analysis (PCA) of Flesh


图9 7种柑橘的类胡萝卜素含量比较  
Comparison of Carotenoid Content of Seven Citrus Fruits