

石墨炉原子吸收法直接测定油品中金属元素——砷和铅

引言


油品中的金属杂质，对加工或油品品质往往造成有害的影响，如：油品中的铅、砷等会导致催化剂中毒。因此，对油品中的痕量元素铅、砷的测定要求很高。目前，主要方法还是需要对油品进行湿法消化，这易导致油品损失，使结果偏低。

本方法采用甲苯稀释，用 PE AA600/800 原子吸收光谱仪，测定时，加钼和硝酸镁作为基体改进剂，同时在石墨炉干燥和灰化阶段中，加入一步通空气技术，消除积碳干扰，直接测定油品中的铅、砷等杂质元素。

实验结果表明：该方法无论在灵敏度、分析速度、准确度和操作等方面，都可以获得非常满意的结果。

仪器

使用 PerkinElmer AAnalyst600/800 原子吸收光谱仪，配以 AS800 自动进样器和 PerkinElmer 专用石墨管，WinLab32 for AA 软件，特殊通空气程序，消




作者

PerkinElmer, Inc.

原理



样品经进样器进样后，在石墨管中通过干燥、碳化、灰化、原子化等步骤，使样品蒸发离解形成原子蒸气，对来自光源的特征电磁辐射产生吸收，将测得的样品吸光度和标准吸光度进行比较，确定样品中被测元素的含量。

优点：无需消化，可以直接分析原油、石脑油、乙醇汽油等。检出限低ppb级，所需样品量少等。



特点


1. THGA(横向加热)石墨炉
2. 专利的纵向交流zeeman效应背景校正
3. 新型内置直流ASC0M石墨炉电源
4. 专利的TTC(真实温度控制)技术,石墨炉温度准确到 $\pm 10^{\circ}\text{C}$,重现性好
5. 气动式锁紧石墨管,石墨管损坏自动报警
6. 管内外气流分开,原子化阶段内管自动停气,氩气耗量 $<0.7\text{L}/\text{min}$
7. 自动节气停水技术
8. 独特的通空气技术,消除有机物、氯化物等的干扰
9. 专利的BOC(自动基线调零)技术,每次测定读数前,自动进行零点漂移校正,克服由于灯、电路漂移或石墨炉两边石英窗变脏引的零点变化,长时间测定稳定性好



准备

石墨炉分析


- 特征质量
- 灵敏度检查
- 基体改进剂配制
- 进样针调节
- 仪器推荐条件
- STPF



约45°, 斜口朝外, 长约0.8cm

斜口尽量与石墨管内口相切

进样针在直径方向上约进入(6~7)/10



标准和样品制备

- 无机标样配制: 称取一定量水标样, 用乙醇稀释至浓度为 50ng/g
- 有机标样配制: 称取一定量有机标样, 用甲苯稀释至浓度为 50ng/g
- 基体改进剂: 称取一定量的 1%Pd 和 1%Mg(NO₃)₂, 用乙醇稀释至浓度为 0.1%Pd 和 0.06%Mg(NO₃)₂
- 稀释液和清洗液: 分析乙醇汽油样品时用乙醇, 分析原油等样品时用甲苯



油品用甲苯稀释5~10倍分析



加氧汽油、石脑油等直接分析

分析条件和结果

1、油品中铅的分析条件和结果

1.1、铅分析条件

WinLab32 AA Furnace - 推荐条件

元素: Pb (Lead)

设置数据

| | |
|-----------|--------------------------------------|
| 波长 (nm) | 283.3 |
| 小狭缝 (nm) | 0.7 |
| 再吸收 (abs) | 1.50 |
| 温度 (°C) | 高温分解: 850 原子化: 1600 |
| 原子化位置 | Pyro/Platform |
| 化学改进剂 | 0.050 mg NH4H2PO4 + 0.003 mg MgND3I2 |

性能检查

特征质量: 典型的 ... 30 pg/0.0044 A-s

灵敏度检查: 50 ug/L for 0.15 A-s

% Modifier = (? mg)(100) / (5 µL)

For the analysis of As with Pd:

% Modifier = (0.005 mg)(100)/(5 µL) = 0.1% Pd

1. 稀释剂: 0.2% 硝酸
2. 由于这一元素巨大的污染问题可能发生
3. 无电极放电管灯能用汞获得这些数据
4. 使用217.0 nm线, 在大约0.5 A之上, 添加溶剂可能发生
5. 推荐原子化时间: 3秒
6. 用 50 µg的特征质量, 这 0.005 mg Pd + 0.003 mg MgND3I2 的修饰剂允许1000 °C的溶解温度, 用这一修饰剂, 更低的1900-°C的氧化温度是必需的
7. 在上述的条件之下, 可能会导致特征质量有 +/- 20%的变化
8. 当使用带堵头的管子时, 特征质量是 21 µg/0.0044 A-s.

方法编辑器: Pb

石墨炉程序




| 步骤 # | 温度 (°C) | 提升时间 | 保持时间 | 内气流 | 气体类型 |
|------|---------|------|------|-----|------|
| 1 | 110 | 5 | 30 | 250 | 标准 |
| 2 | 130 | 15 | 30 | 250 | 标准 |
| 3 | 450 | 15 | 15 | 250 | 指定 |
| 4 | 1000 | 10 | 20 | 250 | 标准 |
| 5 | 1800 | 0 | 5 | 0 | 标准 |
| 6 | 2450 | 1 | 3 | 250 | 标准 |
| 7 | | | | | |
| 8 | | | | | |

读数步骤: 5

进样温度 (°C): 20

额外的石墨炉清洗: 否 是

通空气铅分析条件



Sample Information Editor : Untitled

Parameters Common to All Samples

| | |
|--------------|----|
| Batch ID | |
| Volume Units | mL |
| Weight Units | g |

File Description
Default Sample Information File

*Parameters marked with an asterisk override settings in the method.

Parameters That Vary By Sample

| Sample No | A/S Location | Sample ID | Initial Sample Wt. | Sample Prep. Vol. | Aliquot Volume | Diluted To Vol |
|-----------|--------------|-----------|--------------------|-------------------|----------------|----------------|
| 1 | | | | | | |
| 2 | 12 | 原油 | | | | |
| 3 | 13 | 加氧汽油 | | | | |
| 4 | | | | | | |
| 5 | | | | | | |
| 6 | | | | | | |
| 7 | | | | | | |

Automated Analysis Control

Sample Progress Analysis Progress

Step: Idle

Analyze All Calibrate Analyze Samples Reset Sequence

| Seq. | Loc. | Sample ID | Status |
|------|------|---------------|-------------------|
| 1 | 81 | Calib Blank 1 | Applied |
| 2 | 2 | Calib Std 1 | Applied |
| 3 | 2 | Calib Std 2 | Applied |
| 4 | 2 | Calib Std 3 | Applied |
| 5 | 2 | Calib Std 4 | Applied |
| 6 | 12 | 原油 | Analyzed |
| 7 | 12 | 原油 | Auto Recov Passed |
| 8 | 13 | 加氧汽油 | Analyzed |
| 9 | 13 | 加氧汽油 | Auto Recov Passed |
| 10 | | | |


Method: 1 Rebuild List Print List Priority... Select Loc


Set Up Analyze

1.2 铅结果

1.2.1 乙醇汽油铅图谱和结果

1.2.1.1 乙醇汽油铅图谱





1.2.1.2 乙醇汽油铅结果

```

Results
Sequence No.: 1
Sample ID: Calib Blank 1
Logged In Analyst (Original):
Initial Sample Wt:
Dilution:
Autosampler Location: 81
Date Collected: 2006-11-5 11:55:03
Data Type: Reprocessed on 2006-11-7 22:54:09
Initial Sample Vol:
Sample Prep Vol:

-----
Replicate Data: Calib Blank 1
Repl SampleConc StndConc BlnkCorr Peak Peak Bkgnd Bkgnd Time Peak
# ng/g ng/g Signal Area Height Area ng/g Height
1 [0.00] 0.0043 0.0043 0.0050 0.0051 0.0044 11:58:07 No Stored
Auto-zero performed.

-----
Sequence No.: 2
Sample ID: Calib Std 1
Logged In Analyst (Original):
Initial Sample Wt:
Dilution:
Autosampler Location: 2
Date Collected: 2006-11-5 11:58:03
Data Type: Reprocessed on 2006-11-7 22:54:09
Initial Sample Vol:
Sample Prep Vol:

-----
Replicate Data: Calib Std 1
Repl SampleConc StndConc BlnkCorr Peak Peak Bkgnd Bkgnd Time Peak
# ng/g ng/g Signal Area Height Area ng/g Height
1 [12.5] 0.0360 0.0403 0.0400 0.0122 0.0070 11:58:55 No Stored
Standard number 1 applied. [12.5]
Correlation Coef.: 1.000000 Slope: 0.00288 Intercept: 0.00000

-----
Sequence No.: 3
Sample ID: Calib Std 2
Logged In Analyst (Original):
Initial Sample Wt:
Dilution:
Autosampler Location: 2
Date Collected: 2006-11-5 12:00:51
Data Type: Reprocessed on 2006-11-7 22:54:09
Initial Sample Vol:
Sample Prep Vol:

-----
Replicate Data: Calib Std 2
Repl SampleConc StndConc BlnkCorr Peak Peak Bkgnd Bkgnd Time Peak
# ng/g ng/g Signal Area Height Area ng/g Height
1 [25] 0.0706 0.0748 0.0731 0.0177 0.0125 12:01:43 No Stored
Standard number 2 applied. [25]

```

2/3 Results

Sequence No.: 4 Autosampler Location: 2
Sample ID: Calib Std 3 Date Collected: 2006-11-5 12:03:39
Analyst: Data Type: Reprocessed on 2006-11-7 22:54:09
Logged In Analyst (Original):
Initial Sample Wt: Initial Sample Vol:
Dilution: Sample Prep Vol:

Replicate Data: Calib Std 3

| Repl # | Sample Conc ng/g | Std Conc ng/g | Blk Corr Signal | Peak Area | Peak Height | Bkgnd Area | Bkgnd Height | Time | Peak Stored |
|--------|------------------|---------------|-----------------|-----------|-------------|------------|--------------|------|-------------|
| 1 | [37.5] | 0.1035 | 0.1078 | 0.1022 | 0.0231 | 0.0168 | 12:04:32 | No | |

Standard number 3 applied. [37.5]
Correlation Coef.: 1.000000 Slope: 0.00294 Intercept: 0.00000

Sequence No.: 5 Autosampler Location: 2
Sample ID: Calib Std 4 Date Collected: 2006-11-5 12:06:28
Analyst: Data Type: Reprocessed on 2006-11-7 22:54:09
Logged In Analyst (Original):
Initial Sample Wt: Initial Sample Vol:
Dilution: Sample Prep Vol:

Replicate Data: Calib Std 4

| Repl # | Sample Conc ng/g | Std Conc ng/g | Blk Corr Signal | Peak Area | Peak Height | Bkgnd Area | Bkgnd Height | Time | Peak Stored |
|--------|------------------|---------------|-----------------|-----------|-------------|------------|--------------|------|-------------|
| 1 | [50] | 0.1383 | 0.1426 | 0.1381 | 0.0285 | 0.0229 | 12:07:20 | No | |

Standard number 4 applied. [50]
Correlation Coef.: 0.999823 Slope: 0.00292 Intercept: 0.00000
S-shaped calibration curve detected. Two-coefficient equation used.

Calibration data for Pb 283.3 Equation: Nonlinear Through Zero

| ID | Mean Signal (Abs) | Entered Conc. ng/g | Calculated Conc. ng/g | Standard Deviation | %RSD |
|---------------|-------------------|--------------------|-----------------------|--------------------|------|
| Calib Blank 1 | 0.0000 | 0 | 0.000 | ---- | ---- |
| Calib Std 1 | 0.0360 | 12.5 | 12.555 | ---- | ---- |
| Calib Std 2 | 0.0706 | 25.0 | 24.945 | ---- | ---- |
| Calib Std 3 | 0.1035 | 37.5 | 37.133 | ---- | ---- |
| Calib Std 4 | 0.1383 | 50.0 | 50.369 | ---- | ---- |

Correlation Coef.: 0.999823 Slope: 0.00292 Intercept: 0.00000

2/3 Results

Sequence No.: 6 Autosampler Location: 11
Sample ID: 乙醇汽油 Date Collected: 2006-11-5 12:10:15
Analyst: Data Type: Reprocessed on 2006-11-7 22:54:09
Logged In Analyst (Original):
Initial Sample Wt: Initial Sample Vol:
Dilution: Sample Prep Vol:

Replicate Data: 乙醇汽油

| Repl # | Sample Conc ng/g | Std Conc ng/g | Blk Corr Signal | Peak Area | Peak Height | Bkgnd Area | Bkgnd Height | Time | Peak Stored |
|--------|------------------|---------------|-----------------|-----------|-------------|------------|--------------|----------|-------------|
| 1 | 12.0 | 12.0 | 0.0345 | 0.0387 | 0.0382 | 0.0125 | 0.0075 | 12:11:14 | No |
| 2 | 12.0 | 12.0 | 0.0345 | 0.0388 | 0.0379 | 0.0124 | 0.0091 | 12:14:03 | No |

Mean: 12.0 SD: 0.009 %RSD: 0.075

Sequence No.: 7 Autosampler Location: 11
Sample ID: 乙醇汽油 Date Collected: 2006-11-5 12:15:59
Analyst: Data Type: Reprocessed on 2006-11-7 22:54:10
Logged In Analyst (Original):
Initial Sample Wt: Initial Sample Vol:
Dilution: Sample Prep Vol:

Replicate Data: 乙醇汽油

| Repl # | Sample Conc ng/g | Std Conc ng/g | Blk Corr Signal | Peak Area | Peak Height | Bkgnd Area | Bkgnd Height | Time | Peak Stored |
|--------|------------------|---------------|-----------------|-----------|-------------|------------|--------------|----------|-------------|
| 1 | 24.2 | 24.2 | 0.0686 | 0.0729 | 0.0719 | 0.0185 | 0.0124 | 12:17:01 | No |
| 2 | 24.7 | 24.7 | 0.0700 | 0.0743 | 0.0724 | 0.0177 | 0.0121 | 12:20:01 | No |


Mean: 24.5 SD: 0.357 %RSD: 1.46


Automatic Recovery Check: 乙醇汽油

| Analyte | Expected Conc. | Expected Increase | Measured Increase | Units | Recovery (%) | Limits (%) Lower/Upper |
|----------|----------------|-------------------|-------------------|-------|--------------|------------------------|
| Pb 283.3 | 24.5 | 12.5 | 12.5 | ng/g | 99.9 | 90/110 |

1.2.2 原油铅图谱和结果

1.2.2.1 原油铅图谱





1. 2. 2. 2 原油铅图谱

2/3 Results

Replicate Data: Calib Blank 1

| Repl # | Sample Conc ug/g | Std Conc ug/g | Blk Corr Signal | Peak Area | Peak Height | Bkgnd Area | Bkgnd Height | Time | Peak Stored |
|--------|------------------|---------------|-----------------|-----------|-------------|------------|--------------|------|-------------|
| 1 | [0.00] | 0.0035 | 0.0035 | 0.0042 | 0.0074 | 0.0071 | 13:24:38 | No | |

Auto-zero performed.

Sequence No.: 2
 Sample ID: Calib Std 1
 Analyst:
 Logged In Analyst (Original) :
 Initial Sample Wt:
 Dilution:

Autosampler Location: 2
 Date Collected: 2006-11-5 13:26:35
 Data Type: Reprocessed on 2006-11-7 22:58:23

Initial Sample Vol:
 Sample Prep Vol:

Replicate Data: Calib Std 1

| Repl # | Sample Conc ug/g | Std Conc ug/g | Blk Corr Signal | Peak Area | Peak Height | Bkgnd Area | Bkgnd Height | Time | Peak Stored |
|--------|------------------|---------------|-----------------|-----------|-------------|------------|--------------|------|-------------|
| 1 | [12.5] | 0.0417 | 0.0452 | 0.0396 | 0.0134 | 0.0076 | 13:27:27 | No | |

Standard number 1 applied. [12.5]
 Correlation Coef.: 1.000000 Slope: 0.00334 Intercept: 0.00000

Sequence No.: 3
 Sample ID: Calib Std 2
 Analyst:
 Logged In Analyst (Original) :
 Initial Sample Wt:
 Dilution:

Autosampler Location: 2
 Date Collected: 2006-11-5 13:29:24
 Data Type: Reprocessed on 2006-11-7 22:58:23

Initial Sample Vol:
 Sample Prep Vol:

Replicate Data: Calib Std 2

| Repl # | Sample Conc ug/g | Std Conc ug/g | Blk Corr Signal | Peak Area | Peak Height | Bkgnd Area | Bkgnd Height | Time | Peak Stored |
|--------|------------------|---------------|-----------------|-----------|-------------|------------|--------------|------|-------------|
| 1 | [25] | 0.0816 | 0.0850 | 0.0835 | 0.0222 | 0.0141 | 13:30:16 | No | |

Standard number 2 applied. [25]
 Correlation Coef.: 1.000000 Slope: 0.00342 Intercept: 0.00000

Results

Sequence No.: 4
 Sample ID: Calib Std 3
 Analyst:
 Logged In Analyst (Original):
 Initial Sample Wt:
 Dilution:

Autosampler Location: 2
 Date Collected: 2006-11-5 13:32:13
 Data Type: Reprocessed on 2006-11-7 22:58:24

Initial Sample Vol:
 Sample Prep Vol:

Replicate Data: Calib Std 3

| Repl # | Sample Conc ug/g | Std Conc ug/g | Blk Corr Signal | Peak Area | Peak Height | Bkgnd Area | Bkgnd Height | Time | Peak Stored |
|--------|------------------|---------------|-----------------|-----------|-------------|------------|--------------|------|-------------|
| 1 | [37.5] | 0.1196 | 0.1231 | 0.1229 | 0.0260 | 0.0204 | 13:33:05 | No | |

Standard number 3 applied. [37.5]
 Correlation Coef.: 1.000000 Slope: 0.00342 Intercept: 0.00000

Sequence No.: 5
 Sample ID: Calib Std 4
 Analyst:
 Logged In Analyst (Original):
 Initial Sample Wt:
 Dilution:

Autosampler Location: 2
 Date Collected: 2006-11-5 13:35:02
 Data Type: Reprocessed on 2006-11-7 22:58:24

Initial Sample Vol:
 Sample Prep Vol:

Replicate Data: Calib Std 4

| Repl # | Sample Conc ug/g | Std Conc ug/g | Blk Corr Signal | Peak Area | Peak Height | Bkgnd Area | Bkgnd Height | Time | Peak Stored |
|--------|------------------|---------------|-----------------|-----------|-------------|------------|--------------|------|-------------|
| 1 | [50] | 0.1613 | 0.1648 | 0.1614 | 0.0336 | 0.0272 | 13:35:54 | No | |

Standard number 4 applied. [50]
 Correlation Coef.: 0.999624 Slope: 0.00336 Intercept: 0.00000
 S-shaped calibration curve detected. Two-coefficient equation used.

Calibration data for Pb 283.3 Equation: Nonlinear Through Zero

| ID | Mean Signal (Abs) | Entered Conc. ug/g | Calculated Conc. ug/g | Standard Deviation | %RSD |
|---------------|-------------------|--------------------|-----------------------|--------------------|------|
| Calib Blank 1 | 0.0000 | 0 | 0.000 | ---- | ---- |
| Calib Std 1 | 0.0417 | 12.5 | 12.584 | ---- | ---- |
| Calib Std 2 | 0.0816 | 25.0 | 24.911 | ---- | ---- |
| Calib Std 3 | 0.1196 | 37.5 | 36.966 | ---- | ---- |
| Calib Std 4 | 0.1613 | 50.0 | 50.537 | ---- | ---- |

Correlation Coef.: 0.999624 Slope: 0.00336 Intercept: 0.00000

Results

Sequence No.: 6
 Sample ID: 原油
 Analyst:
 Logged In Analyst (Original):
 Initial Sample Wt:
 Dilution:

Autosampler Location: 12
 Date Collected: 2006-11-5 13:49:36
 Data Type: Reprocessed on 2006-11-7 22:58:24

Initial Sample Vol:
 Sample Prep Vol:

Replicate Data: 原油

| Repl # | Sample Conc ug/g | Std Conc ug/g | Blk Corr Signal | Peak Area | Peak Height | Bkgnd Area | Bkgnd Height | Time | Peak Stored |
|--------|------------------|---------------|-----------------|-----------|-------------|------------|--------------|----------|-------------|
| 1 | 7.49 | 7.49 | 0.0250 | 0.0284 | 0.0575 | 0.0684 | 0.0787 | 13:50:36 | No |
| 2 | 7.54 | 7.54 | 0.0251 | 0.0286 | 0.0590 | 0.0588 | 0.0768 | 13:53:28 | No |
| Mean: | 7.51 | 7.51 | 0.0250 | | | | | | |
| SD: | 0.036 | 0.036 | 0.0001 | | | | | | |
| %RSD: | 0.484 | 0.484 | 0.48 | | | | | | |

Sequence No.: 7
 Sample ID: 原油
 Analyst:
 Logged In Analyst (Original):
 Initial Sample Wt:
 Dilution:

Autosampler Location: 12
 Date Collected: 2006-11-5 13:55:25
 Data Type: Reprocessed on 2006-11-7 22:58:24

Initial Sample Vol:
 Sample Prep Vol:

Replicate Data: 原油

| Repl # | Sample Conc ug/g | Std Conc ug/g | Blk Corr Signal | Peak Area | Peak Height | Bkgnd Area | Bkgnd Height | Time | Peak Stored |
|--------|------------------|---------------|-----------------|-----------|-------------|------------|--------------|----------|-------------|
| 1 | 20.6 | 20.6 | 0.0679 | 0.0713 | 0.1542 | 0.0730 | 0.0980 | 13:58:30 | No |
| 2 | 19.6 | 19.6 | 0.0647 | 0.0682 | 0.1528 | 0.0846 | 0.1054 | 13:59:33 | No |
| Mean: | 20.1 | 20.1 | 0.0663 | | | | | | |
| SD: | 0.699 | 0.699 | 0.0023 | | | | | | |
| %RSD: | 3.47 | 3.47 | 3.40 | | | | | | |

Automatic Recovery Check: 原油

| Analyte | Expected Conc. | Expected Increase | Measured Increase | Units | Recovery (%) | Limits (%) Lower/Upper |
|----------|----------------|-------------------|-------------------|-------|--------------|------------------------|
| Pb 283.3 | 20.0 | 12.5 | 12.6 | ug/g | 101.0 | 90/110 |

Results

Sequence No.: 8
 Sample ID: 加氧汽油
 Analyst:
 Logged In Analyst (Original):
 Initial Sample Wt:
 Dilution:

Autosampler Location: 13
 Date Collected: 2006-11-5 14:12:35
 Data Type: Reprocessed on 2006-11-7 22:58:24

Initial Sample Vol:
 Sample Prep Vol:

Replicate Data: 加氧汽油

| Repl # | Sample Conc ug/g | Std Conc ug/g | Blk Corr Signal | Peak Area | Peak Height | Bkgnd Area | Bkgnd Height | Time | Peak Stored |
|--------|------------------|---------------|-----------------|-----------|-------------|------------|--------------|----------|-------------|
| 1 | 1.55 | 1.55 | 0.0052 | 0.0087 | 0.0154 | 0.0394 | 0.0380 | 14:13:33 | No |
| 2 | 1.45 | 1.45 | 0.0049 | 0.0083 | 0.0145 | 0.0380 | 0.0366 | 14:16:24 | No |
| Mean: | 1.50 | 1.50 | 0.0050 | | | | | | |
| SD: | 0.071 | 0.071 | 0.0002 | | | | | | |
| %RSD: | 4.73 | 4.73 | 4.72 | | | | | | |

Sequence No.: 9
 Sample ID: 加氧汽油
 Analyst:
 Logged In Analyst (Original):
 Initial Sample Wt:
 Dilution:

Autosampler Location: 13
 Date Collected: 2006-11-5 14:18:21
 Data Type: Reprocessed on 2006-11-7 22:58:24

Initial Sample Vol:
 Sample Prep Vol:

Replicate Data: 加氧汽油

| Repl # | Sample Conc ug/g | Std Conc ug/g | Blk Corr Signal | Peak Area | Peak Height | Bkgnd Area | Bkgnd Height | Time | Peak Stored |
|--------|------------------|---------------|-----------------|-----------|-------------|------------|--------------|----------|-------------|
| 1 | 13.9 | 13.9 | 0.0459 | 0.0494 | 0.0856 | 0.0507 | 0.0582 | 14:19:23 | No |
| 2 | 14.6 | 14.6 | 0.0482 | 0.0516 | 0.0893 | 0.0490 | 0.0576 | 14:22:22 | No |
| Mean: | 14.2 | 14.2 | 0.0470 | | | | | | |
| SD: | 0.485 | 0.485 | 0.0016 | | | | | | |
| %RSD: | 3.41 | 3.41 | 3.37 | | | | | | |

Automatic Recovery Check: 加氧汽油

| Analyte | Expected Conc. | Expected Increase | Measured Increase | Units | Recovery (%) | Limits (%) Lower/Upper |
|----------|----------------|-------------------|-------------------|-------|--------------|------------------------|
| Pb 283.3 | 14.0 | 12.5 | 12.7 | ug/g | 101.7 | 90/110 |

2、油品中砷的分析条件和结果

2.1、砷分析条件

WinLab 32 AA Furnace

文件(F) 编辑(E) 工具(T) 分析(A) 选项(O) 窗口(W) 帮助(H)

方法: 无标题
 试样信息: 无标题

推荐条件

元素: As [Arsenic]

设置数据

波长 (nm): 193.7
 小狭缝 (nm): 0.7
 再吸收 (abs): 1.30
 温度 (°C): 高温分解 .. 1200 原子化 .. 2000
 原子化位点: Pyro/Platform
 化学改进剂: 0.005 mg Pd + 0.003 mg Mg/NO3

性能检查

特征质量: 典型的 ... 40 pg/0.0044 A-s
 灵敏度检查: 50 ug/L for 0.11 A-s

% Modifier = (? mg)(100) / (5 µL)
 For the analysis of As with Pd:
 % Modifier = (0.005 mg)(100)/(5 µL)
 = 0.1% Pd

显示推荐条件窗口

方法编辑器: As

石墨炉程序

| 步骤 # | 温度 (°C) | 坡升时间 | 保持时间 | 内气流 | 气体类型 |
|------|---------|------|------|-----|------|
| 1 | 110 | 1 | 30 | 250 | 标准 |
| 2 | 130 | 15 | 30 | 250 | 标准 |
| 3 | 450 | 15 | 15 | 250 | 指定 |
| 4 | 1000 | 20 | 20 | 250 | 标准 |
| 5 | 2100 | 0 | 5 | 0 | 标准 |
| 6 | 2450 | 1 | 3 | 250 | 标准 |
| 7 | | | | | |
| 8 | | | | | |

读数步骤: 5
 进样温度 (°C): 20
 额外的石墨炉清洗: 否 是

默认程序 设置...

石墨炉程序
 自动进样器
 顺序

光谱仪 取样器 校准 检查 OK 质控 选项

通空气铅分析条件

方法编辑器: As

石墨炉自动进样器

试样
 体积 (μL): 20 稀释液体积 (μL): 10
 稀释液位置: 01

基体改进剂
 #1 体积 (μL): 5 位置: 1
 添加校准空白与试样
 添加试剂空白与试样
 #2 体积 (μL): 0 位置: 1
 添加校准空白与试样
 添加试剂空白与试样

基体改进剂
 0.1%Pd+0.06%Mg(NO₃)₂

方法编辑器: As

自动进样器与石墨炉序列

| 步骤 | 操作与参数 |
|----|--|
| A | Pipet diluent + modifier1 + spike + sample/std |
| B | Run furnace steps 1 to end |
| C | |
| D | |
| E | |
| F | |
| G | |
| H | |

修改序列 简单... 自定义...
 吸液速度 (%): 80 取样: 100

石墨炉程序
 自动进样器
 顺序

光谱仪 取样器 校准 检查 OK 质控 选项

取样速度

方法编辑器: As

校准方程式

方程式: 非线性通过零点
 最大小数点位数: 3
 最大有效数字位数: 3
 加入法: 系统准备 分析准备

单位
 校准: ng/g
 试样: ng/g

方法编辑器: As



标样浓度

| 名称 | 浓度 | 进样量 (μL) | 稀释液 (μL) |
|--------|-------------|----------|----------|
| 校准空白 | 校准空白 1 | 81 | 20 |
| 校准空白 | 校准空白 1 | 81 | 20 |
| 校准标准 1 | 校准标准 1 12.5 | 3 | 5 |
| 校准标准 2 | 校准标准 2 25 | 3 | 10 |
| 校准标准 3 | 校准标准 3 37.5 | 3 | 15 |
| 校准标准 4 | 校准标准 4 50 | 3 | 20 |
| 校准标准 5 | | | |
| 校准标准 6 | | | |
| 校准标准 7 | | | |
| 校准标准 8 | | | |

计算标样体积...

方程式与单位
 标样浓度
 初始校准曲线
 校准检查
 重新校准


光谱仪 取样器 校准 检查 OK 质控 选项



2.2 砷结果

2.2.1 乙醇汽油砷图谱和结果

2.2.1.1 乙醇汽油砷图谱



2.2.1.2 乙醇汽油砷结果

Results

Replicate Data: Calib Blank 1

| Repl # | Sample Conc ug/g | Std Conc ug/g | Blk Corr Signal | Peak Area | Peak Height | Bkgnd Area | Bkgnd Height | Time | Peak Stored |
|--------|------------------|---------------|-----------------|-----------|-------------|------------|--------------|------|-------------|
| 1 | [0.00] | -0.0013 | -0.0013 | 0.0045 | 0.0152 | 0.0527 | 11:11:15 | No | |

Auto-zero performed.

Sequence No.: 2
 Sample ID: Calib Std 1
 Analyst:
 Logged In Analyst (Original):
 Initial Sample Wt:
 Dilution:

Autosampler Location: 2
 Date Collected: 2006-11-5 11:13:11
 Data Type: Reprocessed on 2006-11-7 23:02:39
 Initial Sample Vol:
 Sample Prep Vol:

Replicate Data: Calib Std 1

| Repl # | Sample Conc ug/g | Std Conc ug/g | Blk Corr Signal | Peak Area | Peak Height | Bkgnd Area | Bkgnd Height | Time | Peak Stored |
|--------|------------------|---------------|-----------------|-----------|-------------|------------|--------------|------|-------------|
| 1 | [12.5] | 0.0363 | 0.0350 | 0.1118 | 0.0072 | 0.0495 | 11:14:03 | No | |

Standard number 1 applied. [12.5]
 Correlation Coef.: 1.000000 Slope: 0.00291 Intercept: 0.00000

Sequence No.: 3
 Sample ID: Calib Std 2
 Analyst:
 Logged In Analyst (Original):
 Initial Sample Wt:
 Dilution:

Autosampler Location: 2
 Date Collected: 2006-11-5 11:16:00
 Data Type: Reprocessed on 2006-11-7 23:02:39
 Initial Sample Vol:
 Sample Prep Vol:

Replicate Data: Calib Std 2

| Repl # | Sample Conc ug/g | Std Conc ug/g | Blk Corr Signal | Peak Area | Peak Height | Bkgnd Area | Bkgnd Height | Time | Peak Stored |
|--------|------------------|---------------|-----------------|-----------|-------------|------------|--------------|------|-------------|
| 1 | [25] | 0.0707 | 0.0693 | 0.2200 | 0.0112 | 0.0546 | 11:16:52 | No | |

Standard number 2 applied. [25]
 Correlation Coef.: 1.000000 Slope: 0.00299 Intercept: 0.00000

Results

Sequence No.: 4
 Sample ID: Calib Std 3
 Analyst:
 Logged In Analyst (Original):
 Initial Sample Wt:
 Dilution:

Autosampler Location: 2
 Date Collected: 2006-11-5 11:18:48
 Data Type: Reprocessed on 2006-11-7 23:02:39
 Initial Sample Vol:
 Sample Prep Vol:

Replicate Data: Calib Std 3

| Repl # | Sample Conc ug/g | Std Conc ug/g | Blk Corr Signal | Peak Area | Peak Height | Bkgnd Area | Bkgnd Height | Time | Peak Stored |
|--------|------------------|---------------|-----------------|-----------|-------------|------------|--------------|------|-------------|
| 1 | [37.5] | 0.1037 | 0.1023 | 0.3276 | 0.0121 | 0.0568 | 11:19:40 | No | |

Standard number 3 applied. [37.5]
 Correlation Coef.: 1.000000 Slope: 0.00301 Intercept: 0.00000

Sequence No.: 5
 Sample ID: Calib Std 4
 Analyst:
 Logged In Analyst (Original):
 Initial Sample Wt:
 Dilution:

Autosampler Location: 2
 Date Collected: 2006-11-5 11:21:36
 Data Type: Reprocessed on 2006-11-7 23:02:40
 Initial Sample Vol:
 Sample Prep Vol:

Replicate Data: Calib Std 4

| Repl # | Sample Conc ug/g | Std Conc ug/g | Blk Corr Signal | Peak Area | Peak Height | Bkgnd Area | Bkgnd Height | Time | Peak Stored |
|--------|------------------|---------------|-----------------|-----------|-------------|------------|--------------|------|-------------|
| 1 | [50] | 0.1378 | 0.1364 | 0.4400 | 0.0132 | 0.0670 | 11:22:28 | No | |

Standard number 4 applied. [50]
 Correlation Coef.: 0.999856 Slope: 0.00295 Intercept: 0.00000

Calibration data for As 193.7

Equation: Nonlinear Through Zero

| ID | Mean Signal (Abs) | Entered Conc. ug/g | Calculated Conc. ug/g | Standard Deviation | XRSD |
|---------------|-------------------|--------------------|-----------------------|--------------------|------|
| Calib Blank 1 | 0.0000 | 0 | 0.000 | ---- | ---- |
| Calib Std 1 | 0.0363 | 12.5 | 12.584 | ---- | ---- |
| Calib Std 2 | 0.0707 | 25.0 | 24.898 | ---- | ---- |
| Calib Std 3 | 0.1037 | 37.5 | 37.196 | ---- | ---- |
| Calib Std 4 | 0.1378 | 50.0 | 50.345 | ---- | ---- |

Correlation Coef.: 0.999856 Slope: 0.00295 Intercept: 0.00000

Results

Sequence No.: 6
 Sample ID: 乙醇汽油
 Analyst:
 Logged In Analyst (Original):
 Initial Sample Wt:
 Dilution:

Autosampler Location: 11
 Date Collected: 2006-11-5 11:35:05
 Data Type: Reprocessed on 2006-11-7 23:02:40

Initial Sample Vol:
 Sample Prep Vol:

Replicate Data: 乙醇汽油

| Repl # | Sample Conc ug/g | Std Conc ug/g | Blak Corr Signal | Peak Area | Peak Height | Bkgnd Area | Bkgnd Height | Time | Peak Stored |
|--------|------------------|---------------|------------------|-----------|-------------|------------|--------------|----------|-------------|
| 1 | 11.9 | 11.9 | 0.0344 | 0.0331 | 0.1119 | 0.0071 | 0.0514 | 11:38:04 | No |
| 2 | 12.0 | 12.0 | 0.0346 | 0.0333 | 0.1101 | 0.0068 | 0.0489 | 11:38:53 | No |
| Mean: | 11.9 | 11.9 | 0.0345 | | | | | | |
| SD: | 0.053 | 0.053 | 0.0002 | | | | | | |
| %RSD: | 0.447 | 0.447 | 0.44 | | | | | | |

Sequence No.: 7
 Sample ID: 乙醇汽油
 Analyst:
 Logged In Analyst (Original):
 Initial Sample Wt:
 Dilution:

Autosampler Location: 11
 Date Collected: 2006-11-5 11:40:50
 Data Type: Reprocessed on 2006-11-7 23:02:40

Initial Sample Vol:
 Sample Prep Vol:

Replicate Data: 乙醇汽油


| Repl # | Sample Conc ug/g | Std Conc ug/g | Blak Corr Signal | Peak Area | Peak Height | Bkgnd Area | Bkgnd Height | Time | Peak Stored |
|--------|------------------|---------------|------------------|-----------|-------------|------------|--------------|----------|-------------|
| 1 | 24.5 | 24.5 | 0.0696 | 0.0683 | 0.2205 | 0.0098 | 0.0543 | 11:41:52 | No |
| 2 | 24.8 | 24.8 | 0.0703 | 0.0689 | 0.2221 | 0.0045 | 0.0495 | 11:44:52 | No |
| Mean: | 24.6 | 24.6 | 0.0700 | | | | | | |
| SD: | 0.170 | 0.170 | 0.0005 | | | | | | |
| %RSD: | 0.688 | 0.688 | 0.66 | | | | | | |







Automatic Recovery Check: 乙醇汽油

| Analyte | Expected Conc. | Expected Increase | Measured Increase | Units | Recovery (%) | Limits (%) Lower/Upper |
|----------|----------------|-------------------|-------------------|-------|--------------|------------------------|
| As 193.7 | 24.4 | 12.5 | 12.7 | ug/g | 101.7 | 90/110 |

2.2.2 原油砷图谱和结果

2.2.2.1 原油砷图谱





2.2.2.2 原油砷结果

Results

Sequence No.: 1
 Sample ID: 校准空白 1
 Analyst:
 Logged In Analyst (Original):
 Initial Sample Wt:
 Dilution:

Autosampler Location: 81
 Date Collected: 2006-8-25 12:24:23
 Data Type: Reprocessed on 2006-8-26 0:37:40

Initial Sample Vol:
 Sample Prep Vol:

Replicate Data: 校准空白 1

| Repl # | Sample Conc ng/g | Std Conc ng/g | Blk Corr Signal | Peak Area | Peak Height | Bkgnd Area | Bkgnd Height | Time | Peak Stored |
|--------|------------------|---------------|-----------------|-----------|-------------|------------|--------------|------|-------------|
| 1 | [0.00] | [0.00] | -0.0004 | 0.0038 | 0.0269 | 0.0192 | 12:25:17 | No | |
| 2 | [0.00] | 0.0001 | 0.0001 | 0.0026 | 0.0268 | 0.0169 | 12:28:53 | No | |
| Mean: | [0.00] | -0.0001 | | | | | | | |
| SD: | 0.00 | 0.0004 | | | | | | | |
| %RSD: | 0.00 | 270.60 | | | | | | | |

Auto-zero performed.

Sequence No.: 2
 Sample ID: 校准标样 1
 Analyst:
 Logged In Analyst (Original):
 Initial Sample Wt:
 Dilution:

Autosampler Location: 3
 Date Collected: 2006-8-25 12:31:27
 Data Type: Reprocessed on 2006-8-26 0:37:40

Initial Sample Vol:
 Sample Prep Vol:

Replicate Data: 校准标样 1

| Repl # | Sample Conc ng/g | Std Conc ng/g | Blk Corr Signal | Peak Area | Peak Height | Bkgnd Area | Bkgnd Height | Time | Peak Stored |
|--------|------------------|---------------|-----------------|-----------|-------------|------------|--------------|------|-------------|
| 1 | [12.5] | 0.0233 | 0.0232 | 0.0306 | 0.0333 | 0.0165 | 12:32:24 | No | |
| 2 | [12.5] | 0.0241 | 0.0240 | 0.0266 | 0.0268 | 0.0165 | 12:35:53 | No | |
| Mean: | [12.5] | 0.0237 | | | | | | | |
| SD: | 0.0 | 0.0006 | | | | | | | |
| %RSD: | 0.0 | 2.42 | | | | | | | |

Standard number 1 applied. [12.5]
 Correlation Coef.: 1.000000 Slope: 0.00190 Intercept: 0.00000

Results

Sequence No.: 3
 Sample ID: 校准标样 2
 Analyst:
 Logged In Analyst (Original):
 Initial Sample Wt:
 Dilution:

Autosampler Location: 3
 Date Collected: 2006-8-25 12:38:26
 Data Type: Reprocessed on 2006-8-26 0:37:40

Initial Sample Vol:
 Sample Prep Vol:

Replicate Data: 校准标样 2

| Repl # | Sample Conc ng/g | Std Conc ng/g | Blk Corr Signal | Peak Area | Peak Height | Bkgnd Area | Bkgnd Height | Time | Peak Stored |
|--------|------------------|---------------|-----------------|-----------|-------------|------------|--------------|------|-------------|
| 1 | [25] | 0.0474 | 0.0472 | 0.0511 | 0.0333 | 0.0196 | 12:39:24 | No | |
| 2 | [25] | 0.0484 | 0.0483 | 0.0549 | 0.0292 | 0.0193 | 12:42:54 | No | |
| Mean: | [25] | 0.0479 | | | | | | | |
| SD: | 0 | 0.0007 | | | | | | | |
| %RSD: | 0 | 1.49 | | | | | | | |

Standard number 2 applied. [25]
 Correlation Coef.: 1.000000 Slope: 0.00188 Intercept: 0.00000

Sequence No.: 4
 Sample ID: 校准标样 3
 Analyst:
 Logged In Analyst (Original):
 Initial Sample Wt:
 Dilution:



Autosampler Location: 3
 Date Collected: 2006-8-25 12:45:28
 Data Type: Reprocessed on 2006-8-26 0:37:41


Initial Sample Vol:
 Sample Prep Vol:

Replicate Data: 校准标样 3

| Repl # | Sample Conc ng/g | Std Conc ng/g | Blk Corr Signal | Peak Area | Peak Height | Bkgnd Area | Bkgnd Height | Time | Peak Stored |
|--------|------------------|---------------|-----------------|-----------|-------------|------------|--------------|------|-------------|
| 1 | [37.5] | 0.0732 | 0.0731 | 0.0790 | 0.0274 | 0.0192 | 12:46:24 | No | |
| 2 | [37.5] | 0.0724 | 0.0722 | 0.0746 | 0.0260 | 0.0198 | 12:49:54 | No | |
| Mean: | [37.5] | 0.0728 | | | | | | | |
| SD: | 0.0 | 0.0006 | | | | | | | |
| %RSD: | 0.0 | 0.86 | | | | | | | |

Standard number 3 applied. [37.5]
 Correlation Coef.: 1.000000 Slope: 0.00188 Intercept: 0.00000





PerkinElmer, Inc.
 大中华区总部
 地址: 上海张江高科园区李冰路67弄4号
 邮编: 201203
 电话: (021) 3876 9510
 传真: (021) 387 91316
 www.perkinelmer.com.cn



要获取全球办事处的完整列表, 请访问 www.perkinelmer.com.cn/ContactUs

©2009 PerkinElmer, Inc. 保留所有权利。PerkinElmer徽标和外观设计是PerkinElmer的注册商标。HyperDSC 和 PerkinElmer 是PerkinElmer及其子公司在美国或其它国家和地区的注册商标。文中提及的其它非PerkinElmer及其子公司所有的其它商标均为其各自所有者的财产。PerkinElmer保留随时更改此文档的权利, 恕不另行通知。对于编辑、图片或排版错误概不承担任何责任。